T7 RNA Polymerase



Certificate of Analysis

Product information

Lot number: Z072-092024

Storage temperature: -20°C after reconstitution with storage buffer

50 mM Tris-HCl, 100 mM NaCl, 20 mM DTT, 1 mM EDTA, 50% glycerol, 0.1% Triton X-

Storage conditions: 100, pH 7.9

Document version: COA_Z072_092024_v1

Specifications

Assay	Description	Threshold values	Z 072- 92024
Purity	Evaluated by SDS-PAGE, where the integrated peak area of the RNA Polymerase is compared to other protein peaks in the same lane.	>95%	99%
Proteinases	Assessed using Pierce™ Fluorescent Protease Assay Kit (ThermoFisher Scientific) following the manufacturer's guidelines.	<50 ng/ml	11.3
RNases	Assessed using the RNAse Alert kit (Integrated DNA Technologies), following the manufacturer's guidelines.	<100 RFU	84.4
Double-stranded endonucleases	2.5 μ l of the T7 RNAP are incubated with 1 μ g of supercoiled plasmid DNA for 0 hours, 4 hours at 37°C and 16 hours at room temperature and analyzed by agarose gel.	No degradation detected.	Pass
Double-stranded exonucleases	2.5 μ l of the T7 RNAP are incubated with 0.5 μ g of linearized double-stranded DNA substrate for 0 hours, 4 hours at 37°C and 16 hours at room temperature and analyzed by agarose gel.	No degradation detected.	Pass
Single-stranded exonucleases	$2.5\mu l$ of the T7 RNAP are incubated with 1 μg of ssDNA substrate (calf-thymus DNA) for 0 hours, 4 hours at 37°C and 16 hours at room temperature and analyzed by agarose gel.	No degradation detected.	Pass
Full-length transcripts	Analysis of % of full-length RNA transcripts vs abortive transcripts using a Bioanalyzer Agilent RNA 6000 Nano Kit.	>50%	Pass